


Molecular sensor based monitoring of harmful algae

**Johanna Hessel, Kerstin Oetjen and Dr. Katja Metfies
Young investigator group PLANKTOSENS**

**24.02.2015
Session 085**


1. Introduction

How are harmful algae defined?



1. Introduction

How are harmful algae defined?





Toxin group

Okadaic acid

Saxitoxins

Domoic acid




Syndromes

Diarrheal Shellfish Poisoning


Amnestic Shellfish Poisoning

Paralytic Shellfish Poisoning

1. Introduction




Number of records in OBIS of species that can potentially cause a HAB



© modified after UNESCO, HAN No.49/2014

1. Introduction

Detection and monitoring systems



2. Technology approach

Sampling +
Filtration

Cell lysis

Molecular
quantification



Automatic filtration module




Ultrasound module



Nucleic acid biosensor


2.1 Sampling and filtration



Automatic filtration module for marine microbes (AutoFiM)

- Sampling reservoir (up to 5 L)


2.1 Sampling and filtration



Automatic filtration module for marine microbes (AutoFiM)


- Sampling reservoir (up to 5 L)
- Filter magazine

2.1 Sampling and filtration




Automatic filtration module for marine microbes (AutoFiM)

- Samplng reservoir (up to 5 L)
- Filter magazine
- Filtration unit (vacuum pump)



2.2 Ultrasound unit






Further processing:

1. Filter storage (roundel)
2. Ultrasound unit generating cell extracts + free nucleic acids





2.3 Nucleic acid biosensor

Molecular detection and quantification: Sandwich Hybridisation




2.3 Nucleic acid biosensor

Molecular detection and quantification: Sandwich Hybridisation



© modified after Wollschläger et al. 2014

- Redox-reaction measured as electrochemical signal
- Quantification of target species (calibration curves)




3. Main Experiments and Objectives

➤ „Proof of Principle“


(1) Autonomous workflow: combination of AutoFiM + us-unit and biosensor

(2) Evaluation and optimisation of

- Filtration procedure
- Sensor chip design
- Spotting performance
- Signal:noise-ratio
- Species-specific probes
- Calibration curves
- Detection chemicals
- Ultrasound intensity



4. Application and field deployability



© modified after I. Peeken (AWI)

- (1) Cruise Sampling (Polarstern May/June 2015)
 - Identify abundance of potentially toxic dinoflagellates
 - Do distribution ranges of HAB-species expand northwards?

- (2) On-site field sampling (e.g. Orkney Islands)



Polarstern © awi.de

Many thanks to...




Johanna Hessel
Alfred-Wegener-Institut
Am Handelshafen 12
27570 Bremerhaven
Germany
johanna.hessel@awi.de

www.awi.de

www.enviguard.net

Chipdesign and signal:noise-ratio



Multi(8)-carbonchip

